

Steribata®

SURGICAL GOWNS

Double protection with maximum comfort



In the right place,
when life is at risk

Made
in
Europe

FAPOMED

Contacts
Av. Dr. Ribeiro Magalhães, 791
4610-108 Felgueiras
Portugal

T. +351 255 310 680
F. +351 255 310 689
welcome@fapomed.pt

Production Units | GPS
Felgueiras, PT 41° 21.633 N 8° 11.608 W
Baião, PT 41° 9.517 N 8° 2.883 W
Gosha, UA 50° 36.664 N 26° 40.767 W

Accuro®
Sterikit®
Steribata®




QUADRO
DE REFERÊNCIA
ESTRATÉGICO
NACIONAL




ASSOCIADO DE MEMBER OF
Health Cluster Portugal
Polo de Competitividade da Saúde




Steribata®




Small Surgery



Standard



Reinforced



Urology

SURGICAL GOWNS

Designed and produced in accordance with Directive 93/42/EEC, the range of surgical gowns STERIBATA® adapts to use in different types of surgical procedures, providing an effective barrier on the prevention and / or control of infectious agents transmission.

With varying levels of performance in accordance with EN 13795, surgical gowns STERIBATA® respond to all levels of requirement with the maximum comfort and safety for the patient and medical team.

SOFT,
COMFORTABLE AND
LIQUID REPELLENT
MATERIAL



Gowns that allow greater freedom of movement without compromising safety and quality.

PATENTED TECHNIQUE
OF ULTRASOUND UNION
THAT PROVIDES
INCREASED STRENGTH
AND IMPERMEABILITY.

Small Surgery Gowns

Material	SMS 35 grs/m2, green					
Reference code	Size	Front Length	Circumference of the gown	Highest point	Sleeve without cuff	Chest
GOW.100.00	Standard	106	143	112	74	64

Standard Gowns

Material	SMS 35 grs/m2, blue / SMS 35 grs/m2, green / SPUNLACE 68 grs/m2, blue					
Reference code	Size	Front Length	Circumference of the gown	Highest point	Sleeve without cuff	Chest
GOW.1020 / 7020 / 2020	M	111 / 111 / 111	146 / 146 / 146	120 / 120 / 120	72 / 72 / 72	66 / 66 / 66
GOW.1030 / 7030 / 2030	L	121 / 121 / 121	156 / 156 / 156	130 / 130 / 130	77 / 77 / 76	70 / 70 / 70
GOW.1040 / 7040 / 2040	XL	125 / 125 / 125	157 / 160 / 157	135 / 135 / 135	80 / 80 / 82	73 / 73 / 76
GOW.1050 / 7050 / 2050	XXL	141 / 141 / 141	172 / 172 / 172	152 / 152 / 152	85 / 85 / 85	82 / 82 / 82

Reinforced Gowns

Material	SMS 35 grs/m2, blue / SMS 35 grs/m2, green / SPUNLACE 68 grs/m2, blue					
Reference code	Size	Front Length	Circumference of the gown	Highest point	Sleeve without cuff	Chest
GOW.1120 / 7120 / 2120	M	111 / 111 / 111	146 / 146 / 146	120 / 120 / 120	72 / 72 / 72	66 / 66 / 66
GOW.1130 / 7130 / 2130	L	121 / 121 / 121	156 / 156 / 156	130 / 130 / 130	77 / 77 / 76	70 / 70 / 70
GOW.1140 / 7140 / 2140	XL	125 / 125 / 125	157 / 160 / 157	135 / 135 / 135	80 / 80 / 82	73 / 73 / 76
GOW.1150 / 7150 / 2150	XXL	141 / 141 / 141	172 / 172 / 172	152 / 152 / 152	85 / 85 / 85	82 / 82 / 82

Urology Gowns

Material	SMS 35 grs/m2, blue					
Reference code	Size	Front length	Circumference of the gown	Highest point	Sleeve Without cuff	Chest
GOW.3230	L	132	148	142	76	75
GOW.3240	XL	144	146	154	91	73

Wrapped in crepe paper 45gr, high strength, with two towels

Individual Multivac Packaging

Traceable labels in individual packaging

Sterilization by EO

Product complies with the requirements of Directive 93/42 / EEC and its revision 2007/47 / CEE and in accordance with the following standards:

EN ISO 9001; EN ISO 13485;
EN ISO 14971; EN 13795;
ISO 11607; EN 868-5;
EN ISO11737; EN ISO 11135;
EN 1041; EN 980

In the right place,
when life is at risk